



CC152 mouse mAb

Catalog No	YP-mAb-10047
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	CCDC152 Chr5_400
Immunogen	Synthesized peptide derived from human CC152 AA range: 115-165
Specificity	This antibody detects endogenous levels of human CC152
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CC152;CCDC152 Chr5_400;Coiled-coil domain-containing protein 152
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	Detected in stomach.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images