



ELAV3 mouse mAb

Catalog No	YP-mAb-10048
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ELAVL3 HUC PLE21
Immunogen	Synthesized peptide derived from human ELAV3 AA range: 110-160
Specificity	This antibody detects endogenous levels of human ELAV3
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ELAV (embryonic lethal abnormal vision Drosophila) like 3 (Hu antigen C);ELAV like protein 3;ELAV-like protein 3;ELAV3_HUMAN;Elavl3;Embryonic lethal; abnormal vision;drosophila;homolog-like;Hu antigen C;Hu-antigen C;HuC;HUCL; MGC20653;Paraneoplastic cerebellar degeneration associated antigen; Paraneoplastic cerebellar degeneration-associated antigen;Paraneoplastic limbic encephalitis antigen 21;PLE21;ELAV3;ELAVL3 HUC PLE21;ELAV-like protein 3 (Hu-antigen C) (HuC) (Paraneoplastic cerebellar degeneration-associated antigen) (Paraneoplastic limbic encephalitis antigen 21)
Observed Band	40kD
Calculated Molecular Weight	40kD
Cell Pathway	
Tissue Specificity	Brain specific.
Function	
Background	

**matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images