



LDB2 mouse mAb

Catalog No	YP-mAb-10049
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	LDB2 CLIM1
Immunogen	Synthesized peptide derived from human LDB2 AA range: 268-318
Specificity	This antibody detects endogenous levels of human LDB2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CLIM 1;CLIM1;LDB 2;LDB1;LDB2;LIM domain binding 2;LIM domain binding protein 2;LDB2 CLIM1;LIM domain-binding protein 2 (LDB-2) (Carboxyl-terminal LIM domain-binding protein 1) (CLIM-1) (LIM domain-binding factor CLIM1)
Observed Band	
Calculated Molecular Weight	42 kDa
Cell Pathway	Nucleus .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images