



PACS2 mouse mAb

Catalog No	YP-mAb-10050
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	PACS2 KIAA0602 PACS1L
Immunogen	Synthesized peptide derived from human PACS2 AA range: 807-857
Specificity	This antibody detects endogenous levels of human PACS2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	KIAA0602;PACS 2;PACS1L;PACS1-like protein;PACS-2;PACS2;PACS2 KIAA0602 PACS1L;Phosphofurin acidic cluster sorting protein 2 (PACS-2) (PACS1-like protein)
Observed Band	100-130 kDa
Calculated Molecular Weight	98kd
Cell Pathway	Endoplasmic reticulum . Mitochondrion .
Tissue Specificity	Broadly expressed, with greatest levels in skeletal muscle followed by heart, brain, pancreas and testis.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images