



KCNRG mouse mAb

Catalog No	YP-mAb-10054
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	KCNRG CLLD4
Immunogen	Synthesized peptide derived from human KCNRG AA range: 50-100
Specificity	This antibody detects endogenous levels of human KCNRG
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CLLD4;DLTET;Potassium channel regulator;Potassium channel regulatory protein;Protein CLLD4;KCNRG;KCNRG CLLD4;Potassium channel regulatory protein (Potassium channel regulator) (Protein CLLD4)
Observed Band	31 kDa
Calculated Molecular Weight	272 aa, 31 kDa
Cell Pathway	Endoplasmic reticulum .
Tissue Specificity	Ubiquitous in normal tissues and expressed in some tumor tissues.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images