



TRI17 mouse mAb

Catalog No	YP-mAb-10055
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	TRIM17 RBCC RNF16 TERF
Immunogen	Synthesized peptide derived from human TRI17 AA range: 95-145
Specificity	This antibody detects endogenous levels of human TRI17
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	E3 ubiquitin-protein ligase TRIM17;RBCC;RING finger protein 16;RING finger protein terf;RNF16;TERF;Testis RING finger protein;TRI17_HUMAN;TRIM 17; Trim17;Tripartite motif containing 17;Tripartite motif containing protein 17; Tripartite motif protein 17;Tripartite motif-containing protein 17;TRI17;TRIM17 RBCC RNF16 TERF;E3 ubiquitin-protein ligase TRIM17 (EC 6.3.2.-) (RING finger protein 16) (Testis RING finger protein) (Tripartite motif-containing protein 17)
Observed Band	54kD
Calculated Molecular Weight	54kD
Cell Pathway	Cytoplasm . Lysosome .
Tissue Specificity	Almost exclusively in the testis.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images