



VPS28 mouse mAb

Catalog No	YP-mAb-10056
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	VPS28
Immunogen	Synthesized peptide derived from human VPS28 AA range: 19-69
Specificity	This antibody detects endogenous levels of human VPS28
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ESCRT-I complex subunit VPS28;H-Vps28;Vacuolar protein sorting-associated protein 28 homolog;VPS28;Vacuolar protein sorting-associated protein 28 homolog (H-Vps28) (ESCRT-I complex subunit VPS28)
Observed Band	28-30 kDa
Calculated Molecular Weight	25 kDa
Cell Pathway	Cell membrane . Late endosome membrane ; Peripheral membrane protein .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images