



ASIC4 mouse mAb

Catalog No	YP-mAb-10057
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ASIC4 ACCN4
Immunogen	Synthesized peptide derived from human ASIC4 AA range: 96-146
Specificity	This antibody detects endogenous levels of human ASIC4
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ACCN4;Acid sensing ion channel 4;ASIC4;BNAC4;ASIC4 ACCN4;Acid-sensing ion channel 4 (ASIC4) (Amiloride-sensitive cation channel 4) (Amiloride-sensitive cation channel 4;pituitary); pituitary)
Observed Band	70 kDa
Calculated Molecular Weight	647 aa, 70 kDa
Cell Pathway	Membrane ; Multi-pass membrane protein .
Tissue Specificity	Expressed in pituitary gland. Weakly expressed in brain, vestibular system and organ of Corti.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images