



EIF2D mouse mAb

Catalog No	YP-mAb-10061
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	EIF2D HCA56 LGTN
Immunogen	Synthesized peptide derived from human EIF2D AA range: 247-297
Specificity	This antibody detects endogenous levels of human EIF2D
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ligatin;Hepatocellular carcinoma-associated antigen 56;HCA56;EIF2D;EIF2D HCA56 LGTN;Eukaryotic translation initiation factor 2D (eIF2d) (Hepatocellular carcinoma-associated antigen 56) (Ligatin)
Observed Band	70 kDa
Calculated Molecular Weight	584 aa, 65 kDa
Cell Pathway	Cytoplasm .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images