



TPC2A mouse mAb

Catalog No	YP-mAb-10062
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	TRAPPC2 SEDL
Immunogen	Synthesized peptide derived from human TPC2A AA range: 45-95
Specificity	This antibody detects endogenous levels of human TPC2A
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TPC2A;TRAPPC2 SEDL;Trafficking protein particle complex subunit 2 (Sedlin)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm, perinuclear region . Endoplasmic reticulum-Golgi intermediate compartment . Nucleus . Cytoplasm . Localized in perinuclear granular structures. .
Tissue Specificity	Expressed in brain, heart, kidney, liver, lung, pancreas, placenta, skeletal muscle, fetal cartilage, fibroblasts, placenta and lymphocytes.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images