



STAR6 mouse mAb

Catalog No	YP-mAb-10063
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	STARD6
Immunogen	Synthesized peptide derived from human STAR6 AA range: 71-121
Specificity	This antibody detects endogenous levels of human STAR6
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	STAR6;STARD6;StAR-related lipid transfer protein 6 (START domain-containing protein 6) (StARD6)
Observed Band	24kD
Calculated Molecular Weight	24kD
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images