



MPZL3 mouse mAb

Catalog No	YP-mAb-10064
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	MPZL3 UNQ2966/PRO7425
Immunogen	Synthesized peptide derived from human MPZL3 AA range: 59-109
Specificity	This antibody detects endogenous levels of human MPZL3
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	myelin protein zero like 3;Myelin protein zero-like protein 3;UNQ2966/PRO7425; MPZL3;MPZL3 UNQ2966/PRO7425
Observed Band	70 kDa, 28 kDa
Calculated Molecular Weight	235 aa, 26 kDa
Cell Pathway	Membrane ; Single-pass type I membrane protein .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images