



AN36A mouse mAb

Catalog No	YP-mAb-10065
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ANKRD36 ANKRD36A KIAA1641 UNQ2430/PRO499
Immunogen	Synthesized peptide derived from human AN36A AA range: 818-868
Specificity	This antibody detects endogenous levels of human AN36A
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	AN36A;ANKRD36 ANKRD36A KIAA1641 UNQ2430/PRO499;Ankyrin repeat domain-containing protein 36A
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images