



KLH31 mouse mAb

Catalog No	YP-mAb-10066
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	KLHL31 BKLHD6 KBTBD1 KLHL
Immunogen	Synthesized peptide derived from human KLH31 AA range: 342-392
Specificity	This antibody detects endogenous levels of human KLH31
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	KLH31;KLHL31 BKLHD6 KBTBD1 KLHL;Kelch-like protein 31 (BTB and kelch domain-containing protein 6) (Kelch repeat and BTB domain-containing protein 1) (Kelch-like protein KLHL)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	Strongly expressed in skeletal muscle and weakly in heart. According to PubMed:15302408, not expressed in other tissues. According to PubMed:18719355, abundantly expressed in both embryonic skeletal and heart tissues.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images