



CTRB1 mouse mAb

Catalog No	YP-mAb-10067
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	CTRB1 CTRB
Immunogen	Synthesized peptide derived from human CTRB1 AA range: 179-229
Specificity	This antibody detects endogenous levels of human CTRB1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	chymotrypsinogen B1;CTRB;CTRB1;CTRB1 CTRB;Chymotrypsinogen B (EC 3.4.21.1) [Cleaved into: Chymotrypsin B chain A;Chymotrypsin B chain B; Chymotrypsin B chain C]
Observed Band	28 kDa
Calculated Molecular Weight	263 aa, 28 kDa
Cell Pathway	Secreted, extracellular space.
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images