



PITX2 mouse mAb

Catalog No	YP-mAb-10068
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	PITX2 ARP1 RGS RIEG RIEG1
Immunogen	Synthesized peptide derived from human PITX2 AA range: 142-192
Specificity	This antibody detects endogenous levels of human PITX2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	PITX2;PITX2 ARP1 RGS RIEG RIEG1;Pituitary homeobox 2 (ALL1-responsive protein ARP1) (Homeobox protein PITX2) (Paired-like homeodomain transcription factor 2) (RIEG bicoid-related homeobox transcription factor) (Solurshin)
Observed Band	33 kDa
Calculated Molecular Weight	35kD
Cell Pathway	Nucleus.
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images