



RBM38 mouse mAb

Catalog No	YP-mAb-10069
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	RBM38 RNPC1 SEB4
Immunogen	Synthesized peptide derived from human RBM38 AA range: 158-208
Specificity	This antibody detects endogenous levels of human RBM38
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	RNA-binding motif protein 38;RNA binding protein 38;RNA binding motif protein 38;HSRNASEB;CLL-associated antigen KW-5;RBM38;RBM38 RNPC1 SEB4; RNA-binding protein 38 (CLL-associated antigen KW-5) (HSRNASEB) (RNA-binding motif protein 38) (RNA-binding region-containing protein 1) (ssDNA-binding protein SEB4)
Observed Band	25-31 kDa
Calculated Molecular Weight	216 aa, 23 kDa
Cell Pathway	Cytoplasm, cytosol . Nucleus .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images