



ITLN2 mouse mAb

Catalog No	YP-mAb-10071
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ITLN2 UNQ2789/PRO7179
Immunogen	Synthesized peptide derived from human ITLN2 AA range: 114-164
Specificity	This antibody detects endogenous levels of human ITLN2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Intelectin-2;intelectin 2;HL2;HL 2;Endothelial lectin HL-2;ITLN2;ITLN2 UNQ2789/PRO7179;Intelectin-2 (Endothelial lectin HL-2)
Observed Band	36 kDa
Calculated Molecular Weight	36 kDa
Cell Pathway	Secreted .
Tissue Specificity	Expressed only in the small intestine.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images