



TPC3L mouse mAb

Catalog No	YP-mAb-10073
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	BET3L
Immunogen	Synthesized peptide derived from human TPC3L AA range: 95-145
Specificity	This antibody detects endogenous levels of human TPC3L
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TPC3L;BET3L;Trafficking protein particle complex subunit 3-like protein (TRAPPC3-like protein) (BET3-like protein)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Golgi apparatus, cis-Golgi network . Endoplasmic reticulum .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images