



PFD3 mouse mAb

Catalog No	YP-mAb-10075
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	VBP1 PFDN3
Immunogen	Synthesized peptide derived from human PFD3 AA range: 71-121
Specificity	This antibody detects endogenous levels of human PFD3
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	von Hippel-Lindau-binding protein 1;VHL-binding protein 1;VBP-1;PFDN3;PFD3;VBP1 PFDN3;Prefoldin subunit 3 (HIBBJ46) (Von Hippel-Lindau-binding protein 1) (VBP-1) (VHL-binding protein 1)
Observed Band	23-25 kDa
Calculated Molecular Weight	23 kDa
Cell Pathway	Cytoplasm. Nucleus. In complex with VHL can translocate to the nucleus.
Tissue Specificity	Ubiquitous.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images