



SMR3B mouse mAb

Catalog No	YP-mAb-10077
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SMR3B PBII PRL3 PROL3
Immunogen	Synthesized peptide derived from human SMR3B AA range: 1-51
Specificity	This antibody detects endogenous levels of human SMR3B
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	SMR3B;SMR3B PBII PRL3 PROL3;Submaxillary gland androgen-regulated protein 3B (Proline-rich peptide P-B) (Proline-rich protein 3) [Cleaved into: Peptide P-A;Peptide D1A]
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Secreted.
Tissue Specificity	Secreted into saliva by submaxillary gland. Not expressed in heart, brain, lung, liver, skeletal muscle, Kidney, pancreas or placenta.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images