



PCY1B mouse mAb

Catalog No	YP-mAb-10078
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	PCYT1B CCTB
Immunogen	Synthesized peptide derived from human PCY1B AA range: 89-139
Specificity	This antibody detects endogenous levels of human PCY1B
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	PCY1B;PCYT1B CCTB;Choline-phosphate cytidyltransferase B (EC 2.7.7.15) (CCT-beta) (CTP:phosphocholine cytidyltransferase B) (CCT B) (CT B) (Phosphorylcholine transferase B)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	[Isoform 1]: Cytoplasm . Endoplasmic reticulum .; [Isoform 2]: Endoplasmic reticulum .
Tissue Specificity	[Isoform 1]: Highly expressed in testis, placenta, brain, ovary, liver and fetal lung. ; [Isoform 2]: Expressed in brain, liver and fetal lung.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images