



SLAP2 mouse mAb

Catalog No	YP-mAb-10081
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SLA2 C20orf156 SLAP2
Immunogen	Synthesized peptide derived from human SLAP2 AA range: 159-209
Specificity	This antibody detects endogenous levels of human SLAP2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Src like adapter 2;SLAP-2;SLAP2;SLAP 2;MARS;SLA2 C20orf156 SLAP2;Src-like-adapter 2 (Modulator of antigen receptor signaling) (MARS) (Src-like adapter protein 2) (SLAP-2)
Observed Band	28-30 kDa
Calculated Molecular Weight	261 aa, 29 kDa
Cell Pathway	Cytoplasm.; [Isoform 1]: Cell membrane. Cytoplasmic vesicle. Localized to the plasma membrane and intracellular vesicles, including late endosomal vesicles.; [Isoform 2]: Cytoplasm . May be cytoplasmic and is not localized to membranes.
Tissue Specificity	Predominantly expressed in immune system, with highest levels in peripheral blood leukocytes. Expressed in spleen, thymus and lymph nodes. Expressed in T-cells as well as in monocytes, and at low level in B-cells. Also detected in placenta, prostate, skin, retina and colon.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images