



ZNF30 mouse mAb

Catalog No	YP-mAb-10083
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ZNF30 KOX28
Immunogen	Synthesized peptide derived from human ZNF30 AA range: 422-472
Specificity	This antibody detects endogenous levels of human ZNF30
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	KOX28;zinc finger protein 30;ZNF30;ZNF30 protein;ZNF30 KOX28;Zinc finger protein 30 (Zinc finger protein KOX28)
Observed Band	65-71 kDa
Calculated Molecular Weight	623 aa, 71 kDa
Cell Pathway	Nucleus .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images