



CB046 mouse mAb

Catalog No	YP-mAb-10085
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	LINC00299 C2orf46 NCRNA00299
Immunogen	Synthesized peptide derived from human CB046 AA range: 36-86
Specificity	This antibody detects endogenous levels of human CB046
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CB046;LINC00299 C2orf46 NCRNA00299;Putative uncharacterized protein encoded by LINC00299
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images