



MYOZ3 mouse mAb

Catalog No	YP-mAb-10087
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	MYOZ3 FRP3
Immunogen	Synthesized peptide derived from human MYOZ3 AA range: 130-180
Specificity	This antibody detects endogenous levels of human MYOZ3
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Calsarcin 3;CS 3;CS3;FATZ related protein 3;FRP3;MYOZ3;myozenin 3; MYOZ3 FRP3;Myozenin-3 (Calsarcin-3) (FATZ-related protein 3)
Observed Band	33 kDa
Calculated Molecular Weight	251 aa, 27 kDa
Cell Pathway	Cytoplasm, myofibril, sarcomere, Z line . Localized at the Z-line of skeletal muscle. .
Tissue Specificity	Expressed specifically in skeletal muscle. Not detected in heart.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images