



PRIC2 mouse mAb

Catalog No	YP-mAb-10088
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	PRICKLE2
Immunogen	Synthesized peptide derived from human PRIC2 AA range: 145-195
Specificity	This antibody detects endogenous levels of human PRIC2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	EPM5;PRIC2;PRIC2_HUMAN;Prickle homolog 2 (Drosophila);Prickle planar cell polarity protein 2;prickle-like 2 (Drosophila);prickle-like 2;Prickle-like protein 2; PRICKLE2
Observed Band	96kD
Calculated Molecular Weight	96kD
Cell Pathway	Nucleus membrane .
Tissue Specificity	Expressed in brain, eye and testis. Additionally in fetal brain, adult cartilage, pancreatic islet, gastric cancer and uterus tumors.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images