



BTBD3 mouse mAb

Catalog No	YP-mAb-10089
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	BTBD3 KIAA0952
Immunogen	Synthesized peptide derived from human BTBD3 AA range: 382-432
Specificity	This antibody detects endogenous levels of human BTBD3
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	BTB (POZ) domain containing 3;BTB domain containing 3;BTB/POZ domain containing protein 3;BTB/POZ domain-containing protein 3;Btbd3;BTBD3 protein; BTBD3_HUMAN;dJ742J24.1;BTBD3 KIAA0952
Observed Band	58kD
Calculated Molecular Weight	58kD
Cell Pathway	Cytoplasm, cytosol . Nucleus . Translocates from the cytosol to the nucleus in response to neuronal activity. .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images