



DJC16 mouse mAb

Catalog No	YP-mAb-10090
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	DNAJC16 KIAA0962
Immunogen	Synthesized peptide derived from human DJC16 AA range: 171-221
Specificity	This antibody detects endogenous levels of human DJC16
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	DJC16;DNAJC16 KIAA0962;DnaJ homolog subfamily C member 16
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Membrane ; Single-pass type IV membrane protein .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images