



HM20A mouse mAb

Catalog No	YP-mAb-10093
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	HMG20A HMGX1 HMGXB1
Immunogen	Synthesized peptide derived from human HM20A AA range: 208-258
Specificity	This antibody detects endogenous levels of human HM20A
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	FLJ10739;High mobility group 20A;High mobility group protein 20A; HM20A_HUMAN;HMG box containing protein 20A;HMG box domain containing 1;HMG box-containing protein 20A;HMG domain protein HMGX1;HMG domain-containing protein 1;HMG domain-containing protein HMGX1;hmg20a;HMGX1; HMGXB1;Hypothetical protein FLJ10739;iBRAF;Inhibitor of BRAF35;HM20A; HMG20A HMGX1 HMGXB1;High mobility group protein 20A (HMG box-containing protein 20A) (HMG domain-containing protein 1) (HMG domain-containing protein HMGX1)
Observed Band	50 kDa
Calculated Molecular Weight	40kD
Cell Pathway	Nucleus .
Tissue Specificity	Ubiquitous.
Function	
Background	

**matters needing
attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images