



PKHG3 mouse mAb

Catalog No	YP-mAb-10094
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	PLEKHG3 KIAA0599
Immunogen	Synthesized peptide derived from human PKHG3 AA range: 557-607
Specificity	This antibody detects endogenous levels of human PKHG3
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	PKHG3;PLEKHG3 KIAA0599;Pleckstrin homology domain-containing family G member 3 (PH domain-containing family G member 3)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm, cytoskeleton . Colocalizes with actin at the leading edge of polarized cells. .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images