



T53G5 mouse mAb

Catalog No	YP-mAb-10095
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	TP53TG5 C20orf10
Immunogen	Synthesized peptide derived from human T53G5 AA range: 85-135
Specificity	This antibody detects endogenous levels of human T53G5
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	T53G5;TP53TG5 C20orf10;TP53-target gene 5 protein (TP53-inducible gene 5 protein)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm . Nucleus . Cell cycle dependent intracellular localization.
Tissue Specificity	Highly expressed in heart, brain and small intestine. Less abundant in skeletal muscle, spleen, prostate, ovary and colon. A smaller transcript is expressed specifically in the testis.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images