



DPTOR mouse mAb

Catalog No	YP-mAb-10097
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	DEPTOR DEPDC6
Immunogen	Synthesized peptide derived from human DPTOR AA range: 87-137
Specificity	This antibody detects endogenous levels of human DPTOR
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	DEP domain containing 6;DEP domain-containing mTOR-interacting protein; DEP domain-containing protein 6;DEP.6;DEPDC6;DEPTOR;DPTOR_HUMAN; DPTOR;DEPTOR DEPDC6;DEP domain-containing mTOR-interacting protein (DEP domain-containing protein 6)
Observed Band	48kDa
Calculated Molecular Weight	46kD
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images