



KRA43 mouse mAb

Catalog No	YP-mAb-10098
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	KRTAP4-3 KAP4.3 KRTAP4.3
Immunogen	Synthesized peptide derived from human KRA43 AA range: 76-126
Specificity	This antibody detects endogenous levels of human KRA43
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	KRA43;KRTAP4-3 KAP4.3 KRTAP4.3;Keratin-associated protein 4-3 (Keratin-associated protein 4.3) (Ultrahigh sulfur keratin-associated protein 4.3)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	Expressed specifically in the middle/uper portions of the hair cortex. Not detected in the hair matrix or cuticle.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images