



ASXL3 mouse mAb

Catalog No	YP-mAb-10104
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ASXL3 KIAA1713
Immunogen	Synthesized peptide derived from human ASXL3 AA range: 1892-1942
Specificity	This antibody detects endogenous levels of human ASXL3
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ASXL3;ASXL3 KIAA1713;Putative Polycomb group protein ASXL3 (Additional sex combs-like protein 3)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus .
Tissue Specificity	Expressed in pancreatic islets, testis, neuroblastoma, head and neck tumor.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images