



CTGE6 mouse mAb

Catalog No	YP-mAb-10106
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	CTAGE6P CTAGE6
Immunogen	Synthesized peptide derived from human CTGE6 AA range: 169-219
Specificity	This antibody detects endogenous levels of human CTGE6
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CTGE6;CTAGE6P CTAGE6;Putative protein cTAGE-6 (cTAGE family member 6)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Membrane ; Single-pass membrane protein .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images