



CP088 mouse mAb

Catalog No	YP-mAb-10107
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	C16orf88 TSG118
Immunogen	Synthesized peptide derived from human CP088 AA range: 141-191
Specificity	This antibody detects endogenous levels of human CP088
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CP088;C16orf88 TSG118;Protein C16orf88 (Testis-specific gene 118 protein)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus, nucleolus .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images