



F75A1 mouse mAb

Catalog No	YP-mAb-10108
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	FAM75A1 C9orf36
Immunogen	Synthesized peptide derived from human F75A1 AA range: 315-365
Specificity	This antibody detects endogenous levels of human F75A1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	F75A1;FAM75A1 C9orf36;Protein FAM75A1
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Membrane ; Single-pass membrane protein .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images