



MAFIP mouse mAb

Catalog No	YP-mAb-10110
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	MAFIP MIP PP5644
Immunogen	Synthesized peptide derived from human MAFIP AA range: 48-98
Specificity	This antibody detects endogenous levels of human MAFIP
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	MAFIP;MAFIP MIP PP5644;MaFF-interacting protein (Tektin-4 like protein PP5644)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm . Nucleus, nucleolus . Interaction with MAFF induces its translocation from the cytoplasm to the nucleolus.
Tissue Specificity	Strongly expressed in brain, kidney and ovary. Moderately expressed in liver, spleen, thymus, prostate, testis, small intestine and colon. Weakly expressed in heart, placenta, lung and leukocytes.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

