



SRCA mouse mAb

Catalog No	YP-mAb-10112
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SRL
Immunogen	Synthesized peptide derived from human SRCA AA range: 544-594
Specificity	This antibody detects endogenous levels of human SRCA
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	SRCA;SRL;Sarcalumenin
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Sarcoplasmic reticulum lumen . Sarcoplasmic reticulum membrane ; Peripheral membrane protein . May associate with the sarcoplasmic reticulum membrane, via a calcium-dependent mechanism. .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images