



DAZ2 mouse mAb

Catalog No	YP-mAb-10113
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	DAZ2
Immunogen	Synthesized peptide derived from human DAZ2 AA range: 435-485
Specificity	This antibody detects endogenous levels of human DAZ2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	DAZ interacting protein 1/2;DAZ interacting protein testis1;DZIP;DZIP2;DZIPT1; Zinc finger DAZ interacting protein 1;Zinc finger protein DZIP1;Zinc finger protein DZIPT1;DAZ2;Deleted in azoospermia protein 2
Observed Band	110 kDa
Calculated Molecular Weight	99kD
Cell Pathway	Cytoplasm . Nucleus . Predominantly cytoplasmic. Nuclear at some stages of spermatozoide development. Localizes both to the nuclei and cytoplasm of spermatozoide differentiation. Nuclear in fetal gonocytes and in spermatogonial nuclei. It then relocates to the cytoplasm during male meiosis.
Tissue Specificity	Testis specific.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images