



K1161 mouse mAb

Catalog No	YP-mAb-10114
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	KIAA1161
Immunogen	Synthesized peptide derived from human K1161 AA range: 199-249
Specificity	This antibody detects endogenous levels of human K1161
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	K1161;KIAA1161;Uncharacterized family 31 glucosidase KIAA1161 (EC 3.2.1.-)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus membrane ; Single-pass type II membrane protein . Endoplasmic reticulum membrane ; Single-pass type II membrane protein . Only a minor fraction is present in the peripheral endoplasmic reticulum. .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images