



CN142 mouse mAb

Catalog No	YP-mAb-10117
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	C14orf142
Immunogen	Synthesized peptide derived from human CN142 AA range: 22-72
Specificity	This antibody detects endogenous levels of human CN142
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CN142;C14orf142;Uncharacterized protein C14orf142
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images