



S38A6 mouse mAb

Catalog No	YP-mAb-10118
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SLC38A6 NAT1 SNAT6
Immunogen	Synthesized peptide derived from human S38A6 AA range: 27-77
Specificity	This antibody detects endogenous levels of human S38A6
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	N system amino acid transporter NAT-1;N-system amino acid transporter 1;Na(+)-coupled neutral amino acid transporter 6;NAT-1;NAT1;Probable sodium-coupled neutral amino acid transporter 6;S38A6_HUMAN;SLC38A6;SNAT1;Solute carrier family 38 member 6;S38A6;SLC38A6 NAT1 SNAT6;Probable sodium-coupled neutral amino acid transporter 6 (N-system amino acid transporter 1) (NAT-1) (Na(+)-coupled neutral amino acid transporter 6) (Solute carrier family 38 member 6)
Observed Band	51kD
Calculated Molecular Weight	51kD
Cell Pathway	Cell membrane ; Multi-pass membrane protein . Colocalizes with synaptotagmins and SNAP25. .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images