



LAP4A mouse mAb

Catalog No	YP-mAb-10124
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	LAPTM4A KIAA0108 LAPTM4 MBNT MTRP UNQ1846/PRO3574
Immunogen	Synthesized peptide derived from human LAP4A AA range: 131-181
Specificity	This antibody detects endogenous levels of human LAP4A
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	LAP4A;LAPTM4A KIAA0108 LAPTM4 MBNT MTRP UNQ1846/PRO3574; Lysosomal-associated transmembrane protein 4A (Golgi 4-transmembrane-spanning transporter MTP)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Endomembrane system ; Multi-pass membrane protein . May reside in an intracellular membrane-bound compartment. .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images