



NARG2 mouse mAb

Catalog No	YP-mAb-10125
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	NARG2 BRCC1 UNQ3101/PRO10100
Immunogen	Synthesized peptide derived from human NARG2 AA range: 474-524
Specificity	This antibody detects endogenous levels of human NARG2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	BRCC1;ICE2;Interactor of little elongator complex ELL subunit 2;Little elongation complex subunit 2;NMDA receptor-regulated protein 2;NARG2; NARG2 BRCC1 UNQ3101/PRO10100
Observed Band	140 kDa
Calculated Molecular Weight	982 aa, 109 kDa
Cell Pathway	Nucleus . Colocalizes with COIL in subnuclear Cajal and histone locus bodies. Translocates in the LEC complex to Cajal and histone locus bodies at snRNA genes in a ICE1-dependent manner. Associates to transcriptionally active chromatin at snRNA genes.
Tissue Specificity	Expressed at low levels in lung and testis.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images