



FIS1 mouse mAb

Catalog No	YP-mAb-10131
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	FIS1 TTC11 CGI-135
Immunogen	Synthesized peptide derived from human FIS1 AA range: 27-77
Specificity	This antibody detects endogenous levels of human FIS1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CGI 135;CGI-135;FIS1 homolog;hFis1;TPR repeat protein 11;FIS1;FIS1 TTC11 CGI-135;Mitochondrial fission 1 protein (FIS1 homolog) (hFis1) (Tetratricopeptide repeat protein 11) (TPR repeat protein 11)
Observed Band	17 kDa
Calculated Molecular Weight	17 kDa
Cell Pathway	Mitochondrion outer membrane; Single-pass membrane protein. Peroxisome membrane; Single-pass membrane protein.
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images