



PNMA2 mouse mAb

Catalog No	YP-mAb-10134
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	PNMA2 KIAA0883 MA2
Immunogen	Synthesized peptide derived from human PNMA2 AA range: 153-203
Specificity	This antibody detects endogenous levels of human PNMA2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	40 kDa neuronal protein;KIAA0883;MA2;MM2;Onconeural antigen Ma2; paraneoplastic antigen MA2;PNMA2;RGAG2;PNMA2 KIAA0883 MA2; Paraneoplastic antigen Ma2 (40 kDa neuronal protein) (Onconeural antigen Ma2) (Paraneoplastic neuronal antigen MM2)
Observed Band	40 kDa
Calculated Molecular Weight	364 aa, 42 kDa
Cell Pathway	Nucleus, nucleolus .
Tissue Specificity	Brain-specific. In some cancer patients, specifically expressed by testicular tumor cells.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images