



ACAD8 mouse mAb

Catalog No	YP-mAb-10136
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ACAD8 ARC42 IBD
Immunogen	Synthesized peptide derived from human ACAD8 AA range: 326-376
Specificity	This antibody detects endogenous levels of human ACAD8
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ACAD 8;ACAD8;ARC42;IBD;ACAD8 ARC42 IBD;Isobutyryl-CoA dehydrogenase;mitochondrial (EC 1.3.99.-) (Activator-recruited cofactor 42 kDa component) (ARC42) (Acyl-CoA dehydrogenase family member 8) (ACAD-8); mitochondrial (EC 1.3.99.-) (Activator-recruited cofactor 42 kDa component) (ARC42) (Acyl-CoA dehydrogenase family member 8) (ACAD-8)
Observed Band	45 kDa
Calculated Molecular Weight	45 kDa
Cell Pathway	Mitochondrion .
Tissue Specificity	Detected at comparable levels in all tissues examined (heart, lung, brain, skeletal muscle, pancreas and placenta). Weakly expressed in liver and kidney.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images