



SNX20 mouse mAb

Catalog No	YP-mAb-10137
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SNX20 SLIC1
Immunogen	Synthesized peptide derived from human SNX20 AA range: 90-140
Specificity	This antibody detects endogenous levels of human SNX20
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	SLIC 1;SLIC1;SNX20;sorting nexin 20;SNX20 SLIC1;Sorting nexin-20 (Selectin ligand-interactor cytoplasmic 1)
Observed Band	36 kDa
Calculated Molecular Weight	316 aa, 36 kDa
Cell Pathway	Early endosome membrane ; Peripheral membrane protein ; Cytoplasmic side . Cell membrane . Cytoplasm . Nucleus .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images